

Beauty white.



Jet Air Polisher

QUICK JET M



Polishing



Autoclavable



Compact

■ Your patients want white teeth.


The powder jet stream readily removes tobacco tar, tea tannin, stain, and plaque. With Quick Jet M, tooth whiteness and brilliance will be recovered speedily.

Compact design, handpiece type.



Easily portable space-saving design for storage. Attachable to turbine hose and no complicated installation is required. The user-friendly designed grip fits more comfortably in your hand.

Available in quick-disconnect couplings of major brands.



Quick Jet M is available various models in quick-disconnect couplings of major handpiece manufacturers as well as in 2-Hole and 4-Hole connectors to meet demands of users.



More cleaning efficiency with newly designed head.

The smaller nozzle outlet enables pinpointed and efficient cleaning in a short time. The rotary head easily cleans areas hard-to-reach with conventional jet polishers.



Sterilizable in an autoclave.



All components, not only the head are sterilizable in an autoclave. You can use it safely and hygienically. Inner powder passages can be cleaned from tip to toe easily and quickly.

Available in 2 colors.



Quick Jet M comes in 2 colors.
 ○ **ivory white**
 ● **black**
 You can choose the preferred and suitable body color for your dental office.

Specifications

| Model | QUICK JET M - K | QUICK JET M - W | QUICK JET M - S | QUICK JET M - B | QUICK JET M - N | QUICK JET M - 4H | QUICK JET M - 2H |
|---------------------------|---|----------------------|------------------------|-------------------|--------------------------|---------------------|---------------------|
| Connection | MULTiflex® coupling | Roto Quick® coupling | Sirona® Quick coupling | Unifix® coupling | Mach®/Phatelus® coupling | ISO standard 4-Hole | ISO standard 2-Hole |
| Powder consumption volume | 4 g/min at coupling air pressure 0.4 MPa | | | | | | |
| Air pressure | 0.25 to 0.4 MPa (36 to 58 psi, 2.5 to 4.0 bar, 2.5 to 4.1 kgf/cm ²) | | | | | | |
| Length | 201 mm (7.9 inch) | 200 mm (7.9 inch) | 196 mm (7.7 inch) | 196 mm (7.7 inch) | 208 mm (8.2 inch) | 207 mm (8.1 inch) | 204 mm (8.0 inch) |
| Weight | 142 g (5.0 oz.) | 136 g (4.8 oz.) | 137 g (4.8 oz.) | 133 g (4.7 oz.) | 144 g (5.1 oz.) | 143 g (5.0 oz.) | 143 g (5.0 oz.) |

MULTiflex® is a registered trademark of KaVo Dental GmbH.
 Roto Quick® is a registered trademark of W&H Deutschland GmbH.
 Sirona® is a registered trademark of Sirona Dental Systems GmbH.

Unifix® is a registered trademark of Bien-Air Holding SA.
 Mach and Phatelus® are registered trademarks of NAKANISHI INC.
 Specifications are subject to change without prior notice for improvements.

Manufacturer : MICRON CORPORATION

2-17-7, Ikegami, Ota-ku, Tokyo, 146-0082, Japan
 Tel : +81 3 3755 0396, Fax : +81 3 5747 5396
<http://www.micdent.com/> E-mail : contact@micdent.com

European Authorized Representative : Advena Ltd.

Tower Business Centre, 2nd Flr., Tower Street, Swatar, BKR 4013 Malta



<https://premier-dental.com.ua>

