

TENEO. TRANSCENDENTAL

TREATMENT CENTER



WHEN YOUR DENTISTRY IS EXTRAORDINARY, YOUR TREATMENT CENTER IS TRANSCENDENTAL.

TENEO is constantly evolving to meet your needs today, tomorrow, and always. This makes it a logical choice for practitioners who want to work now and in the future with the most modern and comfortable technology while maintaining the highest standards of quality, design, and workflow. With a variety of equipment options for a completely customizable operator, TENEO is an investment in the future of your practice.



TRANSCENDENTAL

TENEO, an entirely new generation of treatment centers, eclipses every expectation with respect to design, comfort, and technology. And as the vanguard of digital dentistry, TENEO also represents the epitome of innovation, combining superior esthetics and premium performance that elevates any practice with its presence.

| INTEGRATION



TENEO IS THE HUB OF YOUR PRACTICE.

With freely programmable settings and a selectable EasyMode, the EasyTouch menu guidance can easily be adapted to your usual working methods. Built-in implantology and endodontic functions make tabletop units and long setup times a thing of the past. You can also control SIDEXIS 4, Windows Media Player, PowerPoint, and more directly from the EasyTouch on the delivery unit.

The CEREC Omnicam system is a great CAD/CAM solution that can be integrated with our Transcendental treatment centers. The CEREC AI (Acquisition Integrated) system integrates into the TENEO treatment unit via a separate arm and a PC. Integrating the HELIODENT^{PLUS} on the light post allows you to take intra-oral x-ray images with your Schick sensor. This solution, or the CEREC AI, can be used without leaving the treatment center.



INTEGRATED IMPLANTOLOGY

No need for a separate tabletop unit—with TENEo, implantology is part of the package—allowing you to seamlessly incorporate your GALILEOS-CEREC integration workflow. Using quick-adjustment keys along with hot keys assigned to the functions you use most often, you can set up your workflow to perform treatments in an optimal, time-saving manner.



INTEGRATED ENDODONTICS

The integrated endodontics function makes tabletop units unnecessary and, along with precise torque control, also offers an ApexLocator for the highest degree of clinical safety. An extensive file library is included for most Endo file systems available. Personalization is also possible by saving your own libraries, sequences, and settings.



I DESIGN

PURE, MINIMALISTIC DESIGN FOR OPTIMAL FORM AND FUNCTION

A winner of several prestigious international design awards, TENE0 features a perfect combination of form and function. The motorized track delivery automatically glides to the determined position for convenient and efficient treatment. Clean lines, smooth surfaces, and premium materials make for a simple yet sophisticated appearance that is at once both stylish and practical. The luxurious patient chair with innovative elevation system, lounge upholstery, massage function, and active lumbar support are a feast for the eyes as well as the senses.

Operated by the foot control, the motor-driven sliding track enables you to effortlessly move the dentist element to the correct position. The motorized headrest is automatically adjusted to the saved patient settings*, providing optimal comfort for them, and optimal access for you.

*Available with SIDEXIS 4.





| EXPERIENCE

A PREMIUM TREATMENT CENTER DESIGNED FOR FIRST-CLASS DENTISTRY.

You provide your patients with an incomparably high level of treatment—that is why they come to you. They realize that your practice offers the comfort, convenience, and comprehensive range of services that no other office can match. You've set their expectations, now it's time to exceed them with TENEO.



TENEO PUTS YOUR WORKFLOW IN MOTION.

TENEO features integrated tabletop device functions, intuitive handling, and innovative controls with many additional options that allow you to offer treatment efficiently, quickly, and comfortably. Benefit from incomparable comfort as a practitioner, while your patients are as relaxed as ever.

HIGHEST POSSIBLE COMFORT FOR PROFESSIONAL AND PATIENT

Ergonomically designed with smart integrated functions, TENEO ensures healthy and comfortable work, so you can achieve optimal results from your first to your last patient.

INTUITIVE USER INTERFACE

Easy to use with customizable configurations that complement the way you work.

INTEGRATION OF DENTAL SPECIALTIES

Facilitates efficient endodontics and implantology treatment without the need for separate tabletop units.

TENEO EQUIPMENT OPTIONS

OPTIONS FOR THE PATIENT CHAIR

Patient-specific chair positions*	<input checked="" type="checkbox"/>
Motor-driven headrest with automatic patient specific headrest adjustment*	<input checked="" type="checkbox"/>
Lounge upholstery	<input checked="" type="checkbox"/>
Lumbar support / massage function	<input checked="" type="checkbox"/>
Right armrest can be swiveled aside	<input checked="" type="checkbox"/>
Left armrest	<input checked="" type="checkbox"/>
4 way foot switch	<input checked="" type="checkbox"/>
Advanced electronic foot switch, cable version	<input checked="" type="checkbox"/>
Booster cushion for small patients	<input type="checkbox"/>

WATER UNIT

Cleaning device with water sanitization function	<input checked="" type="checkbox"/>
Suction hose cleaning with chemicals	<input checked="" type="checkbox"/>
Cuspidor	<input type="checkbox"/>
Purge function with autopurge	<input checked="" type="checkbox"/>

Included Optional

*Available with SIDEXIS 4.

** HELIODENT^{PLUS} or CEREC AI cannot be combined.

DENTIST ELEMENT

SPRAYVIT M with light	<input checked="" type="checkbox"/>
BL ISO C electric motor	<input checked="" type="checkbox"/>
BL Implant electric motor	<input checked="" type="checkbox"/>
Additional BL ISO C electric motor	<input type="checkbox"/>
Turbine hose	<input checked="" type="checkbox"/>
SIROSONIC TL ultrasonic handpiece	<input checked="" type="checkbox"/>
Implantology function	<input checked="" type="checkbox"/>
Endodontic treatment function	<input checked="" type="checkbox"/>
ApexLocator with display on SIVISION monitor*	<input checked="" type="checkbox"/>
Integrated pump for saline solution	<input checked="" type="checkbox"/>
USB interface	<input checked="" type="checkbox"/>
Hands-free operation with cursor control	<input checked="" type="checkbox"/>
Contra-angle handpiece T1 Line ENDO 6 L with light	<input checked="" type="checkbox"/>

ASSISTANT ELEMENT

Surgical suction device	<input checked="" type="checkbox"/>
SPRAYVIT M with light	<input checked="" type="checkbox"/>
Suction device	<input checked="" type="checkbox"/>
Saliva ejector	<input checked="" type="checkbox"/>
Suction hose with slide	<input checked="" type="checkbox"/>

SIVISION DIGITAL

SIROCAM AF+	<input checked="" type="checkbox"/>
22" monitor on the light support system	<input checked="" type="checkbox"/>
Cable set for PC connection	<input checked="" type="checkbox"/>
SIDEXIS 4	<input checked="" type="checkbox"/>

STOOLS

1st stool, HUGO or CARL Manual	<input checked="" type="checkbox"/>
2nd stool, HUGO or CARL Manual Plus	<input checked="" type="checkbox"/>

ADDITIONAL OPTIONS

Tray with support arm	<input type="checkbox"/>
PC interface (Ethernet)	<input checked="" type="checkbox"/>
PC with EasyTouch	<input checked="" type="checkbox"/>
Remote diagnosis	<input checked="" type="checkbox"/>
LEDview operating light	<input checked="" type="checkbox"/>
HELIODENT ^{PLUS} mounted on light post**	<input type="checkbox"/>
CEREC AI mounted on light post**	<input type="checkbox"/>





CEREC AI, (Acquisition Integrated), camera with arm, PC and cradle



HUGO MANUAL (without foot ring)



TRAY WITH SUPPORT ARM



HUGO MANUAL PLUS



HELIODENT PLUS



CUSPIDOR



BOOSTER CUSHION



SIROCAM AF+



22" HD, wide screen monitor

SIRONAUSA.COM

CALL YOUR REPRESENTATIVE FOR MORE INFORMATION.

Your Representative will help answer any questions and guide you through the process of getting started with TENE0.

+38 032 259 04 44

+38 097 445 66 73

SIRONA.KIEV.UA

PREMIER-DENTAL
79020, UKRAINE. LVIV
69 Varshavska str., buil."B"



sirona.
The Dental Company